

## SWS62221 Spectral Sensor

### General Description

NeoSpectra™ sensors are the most compact and the lowest cost Fourier Transform InfraRed (FT-IR) spectrometers on the market today. They deliver the same functionality as standard bench top FT-IR spectrometer instruments: quantification and identification of materials by measuring their spectral response. NeoSpectra's core technology is based on semiconductor Micro Electro Mechanical Systems (MEMS) microfabrication techniques promising unprecedented economies of scale.

The SWS62221 sensor determines the spectral content of the input light in NIR range (selectable options between 1,250 – 2,500 nm). The light is coupled via an optical fiber, and spectrum data are transferred to external host computer via USB cable.

### Features

- Low cost embedded NIR spectral sensor solution
- Smallest FT-IR solution with a single photodetector
- Wide spectral ranges ( $\lambda$ : 1,250 – 2,500 nm)
- Designed for high volume production with economies of scale
- Fast, on-chip, data processing
- Alignment free optics
- User selectable resolutions
- Low power consumption
- Competitive performance


### Applications

The SWS62221 is a spectral sensor that is designed to be used as an enabler of application specific spectroscopy solutions in different application areas and industries.

Typical application areas include:

- Food and agriculture
- Pharmaceutical
- Oil and gas
- Polymers
- Biotech
- Metals
- Industrial
- Chemicals

### Block Diagram



## Specifications

Parameter	Conditions	SWS62221-1.7	SWS62221-2.1	SWS62221-2.5	Units
Wavelength Range	PSD <sup>a</sup> > max PSD/10	1,250 - 1,700	1,300 - 2,100	1,350 - 2,500	nm
Resolution	At $\lambda=1,550$ nm, FWHM criterion	8 or 16			nm
		33.3 or 66.6			cm <sup>-1</sup>
Typical SNR (rms)	2 s Scan time, Tx <sup>b</sup> , resolution = 16 nm <sup>c</sup>	>3000:1 (@ $\lambda$ = 1,600 nm)	> 3000:1 (@ $\lambda$ = 1,950 nm)	3000:1 (@ $\lambda$ = 2,050 nm)	-
	2s Scan time, Rx <sup>d</sup> , resolution = 16 nm <sup>e</sup>	>1,000:1 (@ $\lambda$ = 1,600 nm)	>1,000:1 (@ $\lambda$ = 1,950 nm)	1,000:1 (@ $\lambda$ = 2,050 nm)	-
Temperature	Operation	-5: 40 <sup>f</sup>			°C
Wavelength Accuracy	@ $\lambda$ = 1,400 nm; temperature < 40°C	± 1.5			nm
Wavelength Repeatability	@ $\lambda$ = 1,400 nm; absorbance level = 0.5 A.U., Resolution: 16 nm	± 0.1			nm
Power Consumption	USB powered with 5 V supply	750			mW

## Interfaces

Electrical Interface	USB 2.0
Optical Interface	FCPC Multimode optical fiber, Core diameter $\geq$ 400 $\mu$ m, NA = 0.22
Software Interface	Software: SpectroMOST Basic Edition & SDK Edition OS: Windows XP, Vista, 7, 8, Linux Ubuntu 12.04, Linux Debian Machines: 32 & 64 bits

## Typical Setups

	Transmission Setup (TX)	Reflection Setup (RX)
Typical Coupled Power Value within wavelength range	> 3 mW	> 100 $\mu$ W

<sup>a</sup> PSD: Power Spectral Density – Single beam spectrum

<sup>b</sup>Tx: Typical transmission setup



<sup>c</sup> Performance varies linearly with resolution

<sup>d</sup> Rx: Typical reflection setup

<sup>e</sup> Performance varies linearly with resolution

<sup>f</sup> Can be extended to 70°C upon request

### Specifications and parameters definitions

Parameter	Definition
Wavelength range	The wavelength range is defined as the range where the spectral data are useful. The upper and lower wavelength limits are determined by the wavelength points where the power spectral density reaches one tenth of maximum power spectral density over the range.
Typical SNR	SNR is calculated from the root mean square noise ( $N_{rms}$ ), which is the standard deviation of 100 consecutive 100% lines at each wavelength. $SNR = 1/N_{rms}$
Resolution	<p>Resolution is defined as the minimum spacing between two consecutive wavelength (<math>\Delta\lambda</math>) / wavenumber (<math>\Delta\nu</math>) points that can be fully resolved by the module.</p> <p>Two consecutive lines are fully resolved if separation &gt; Full Width Half Maximum (FWHM) power density of either line.</p> <p>The resolution in wavenumber is constant across the spectral range. The relationship between the resolution in wavelength <math>\Delta\lambda</math>, and the resolution in wavenumber <math>\Delta\nu</math> is governed by <math>\Delta\lambda = \Delta\nu \lambda^2</math>.</p> <div style="display: flex; justify-content: space-around;"> <div style="text-align: center;">  <p>(a)</p> </div> <div style="text-align: center;">  <p>(b)</p> </div> </div> <p style="text-align: center;"><b>Typical variation of wavelength resolution across the spectral range</b>  <b>(a) <math>\Delta\lambda = 8 \text{ nm}</math> (b) <math>\Delta\lambda = 16 \text{ nm}</math></b></p>
Wavelength accuracy	Wavelength accuracy is the difference between the measured wavelength of a wavelength standard (e.g. liquid methylene chloride), and the nominal wavelength reported for that wavelength standard.

### Revision History

Revision	Date	Description
1.0	6/22/15	Initial version

### Contact Information

<p><b>Cairo Office (Headquarters)</b> Si-Ware Systems 3, Khaled Ibn Al-Waleed Street Sheraton, Heliopolis Cairo 11361 Egypt Tel.: +20 2 22 68 47 04 Fax: +20 2 22 68 47 05</p>	<p><b>Los Angeles Office</b> Si-Ware Systems, Inc. 1150 Foothill Blvd., Suite M La Canada, CA 91011 USA Tel.: +1 818 790 1151</p>
--	---

Information in this document is provided in connection with Si-Ware Systems products. These materials are provided by Si-Ware Systems as a service to its customers and may be used for informational purposes only. Si-Ware Systems assumes no responsibility for errors or omissions in these materials. Si-Ware Systems may make changes to its products, specifications, and product descriptions at any time, without notice. Si-Ware Systems makes no commitment to update the information and shall have no responsibility whatsoever for conflicts, incompatibilities, or other difficulties arising from future changes to its products and product descriptions. No license, express or implied, by estoppels or otherwise, to any intellectual property rights is granted by this document. Except as may be provided in Si-Ware Systems' Terms and Conditions of Sale for such products, Si-Ware Systems assumes no liability whatsoever.